RESPONDER PELVIC BINDER

INSTRUCTIONS FOR USE


PRECAUTIONS

- Remove objects from patient's pockets and clothing prior to application.
- Monitor the patient's skin for pressure injury while in use.
- Reuse of device may cause illness (e.g., infection) or injury (reduced performance).
- Dispose of used device per healthcare facility biohazard protocol.

1.) Locate one of the tear tabs on any of the four corners of the packaging. Tear package open, all the way across. Remove Pelvic Binder from vacuum sealed storage bag.

2.) Unroll and slide Pelvic Binder under supine patient. Center the Pelvic Binder over greater trochanters. Verify its placement carefully.

3.) Cut or fold back the free end of the binder to a length that allows for a 5-6 in gap between the two edges when brought together.

4.) Lift the velcro plate across the patients pelvis and place on the opposite side. Gently press velcro together to fully engage.

RESPONDER PELVIC BINDER

INSTRUCTIONS FOR USE

5.) Tighten by pulling both straps while pulling pad handle in opposite direction as counter force; until patient is stabilized. Adjust both velcro tabs individually to achieve form fitting pelvic reduction. Wrap the velcro tabs around the binder. Optional: cut excess velcro not being used.

SECURE VELCRO SAFETY FLAP

Once the binder is placed correctly and tension has been applied, secure the velcro safety strap. This will prevent an accidental release of the Binder Pad.

FINAL VIEW

CONSIDERATIONS:

- Hip circumference?- Responder Binders are 5ft in length so they can be cut to fit for each patient.
- Disposable?- Yes
- Usage?- The Responder Binder is typically recommended for temporary pelvic stabilization, 24-48 hrs use max.
- Skin check?- Daily, routine assessments should be performed.

INTENDED USE: The Responder Binder is a non-invasive, circumferential pelvic belt intended to stabilize pelvic fractures during transport to a definitive care facility.

INDICATIONS FOR USE: The Responder Binder is indicated for application by prehospital and medical professionals or trained personnel, to individuals who have, or are suspected of having, a pelvic ring fracture where hemorrhage might occur.

CONTRAINDICATIONS: Do not utilize the Responder Binder if anatomical landmarks are absent.

SPECIFICATIONS:

- The pelvis can be fixed with appropriate tension by pressing the pelvis from the outer circumference at the height of the greater trochanter of the lower buttocks.
- Available in BLACK,CAMO & BLUE.

INSTRUCTIONS FOR USE
RESPONDER PELVIC BINDER
PRECAUTIONS

Remove objects from patient's pockets and clothing prior to application.

Monitor the patient's skin for pressure injury while in use.

Device is single use only.

Reuse of device may cause illness (e.g., infection) or injury (reduced performance).

Dispose of used device per healthcare facility biohazard protocol.

- 1.) Remove Pelvic Binder from storage bag. Pinch and slide tab to release the top.
- 2.) Unroll and slide pelvic binder under supine patient. Center the pelvic binder over greater trochanters.
- 3.) Cut or fold back free end of binder leaving 4-5 inches between the two edges.
- 4.) Pull the velcro plate across the patient's pelvis and place on the opposite side. Gently press to engage velcro attachment.
- 5.) Pull on both velcro tabs, individually or at the same time to achieve pelvic reduction. Wrap the velcro tabs around the binder over their respective lanes. Optional: cut excess velcro not being used.

INTENDED USE: The Responder Binder is a non-invasive, circumferential pelvic belt intended to stabilize pelvic fractures during transport to a definitive care facility.

INDICATIONS FOR USE: The Responder Binder is indicated for application by prehospital and medical professionals or trained personnel, to individuals who have, or are suspected of having, a pelvic ring fracture where hemorrhage might occur.

CONTRAINDICATIONS: Do not utilize the Responder Binder if anatomical landmarks are absent.

SPECIFICATIONS:

The pelvis can be fixed with appropriate tension by pressing the pelvis from the outer circumference at the height of the greater trochanter of the lower buttocks.

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Verfügbar a BLACK, RED, CAMO & BLAU.

SAFE DISPOSAL: The product should be disposed in the dedicated area as recommended by the local authorities.

NOTICE: Please report to FareTec or to the competent authority of the Member State in which you use this product any serious incident that has occurred and you suspect is in relation to this device. If you suspect a potential defect or a product has not met your expectations, please report this to the Company and provide as much information as possible. Company phone number: +1440-350-9510

FareTec Inc, 1610 West Jackson St, Painesville, OH 44077 USA, (p) 440350-9510
(f) 440-350-9520, www.faretec.com

CMC Medical Devices & Drugs S.L., C/ Horacio Lengo n18, C.P 29006
Málaga-Spain

	Manufacturer		Authorized representative in the European Community
	Medical Device		Importer
	Consult instructions for use		Catalogue number
	Date of manufacture		Batch code