

CT-EMS

Carbon Traction Leg Splint



“A Splint to solve your leg traction needs.. in a small one pound bag.”

4:1 Traction Pulley System has revolutionized traction by enabling precise traction application on all sized patients.

Allows for easy Micro-Adjustments



- ▲ Entire Splint, and ankle hitch, are designed to fit Adult through Pediatric sized Patients.
- ▲ Integral Reflective Strap System; comes out of the bag ready to go, no loose straps.
- ▲ Requires minimal patient movement with no Perineal nerve impingement.
- ▲ Tested extensively in the field by numerous organizations through numerous areas of the world.