


CT-6
VS
KTD



Features and Benefits:

4:1 Adjusting System: This produces more power than the KTD strap system. Also it allows for a micro adjustment to generate a precise amount of traction. The easy locking system means that subsequent adjustment during transport is easy to apply. KTD uses a piece of poly webbing and a buckle for adjustment.



Military Black for minimal visibility, carbon and Zytel offer MINIMUM RADAR Reflectivity.


Compact Size - takes up less pack space.



CT-6 = 360 grams KTD = 622 grams

Weight Savings - Modern construction techniques produce a splint that is LIGHTER than the competition while at the same time being stronger and more impact resistant.

Design - The CT-6 is the most advanced traction splint available today. Improvements in design and construction make this splint the correct choice for tough field application where top performance is expected at all times, under all conditions.



CT-6 fits Adult and Pediatric