

CT-6 Carbon Traction Leg Splint



High Tensile Carbon Tubes-Strong And Stiff



New Ischial Strap:

- Uses a forged aluminum buckle.
- 1.5" in Nylon webbing.
- 2" Elastic patient comfort strap.



4:1 Pulley Traction System:

- Makes micro adjustments easy.
- Built in cleat to hold traction.
- Creates sufficient traction on the largest of patients.
- New Safety Clip to prevent accidental release.

- ✦ Ankle Hitch Designed for all sizes.
- ✦ Field Tested: 30,000+ units in use.
- ✦ No Perineal Nerve Impingement.
- ✦ Minimal patient movement.
- ✦ Extender Leg Strap, for large patients.
- ✦ New Carbon tube construction for added durability.



Unit designed to fit Adult and Pediatric Patients.
All units now include Kwikpoint instruction manuals! Training Video available online.

Unit weighs 500grams
Fits in small bag: 27x8x8cm.
NSN# 6515 01 521 5730
Part # 1126624-M

BoundTree
medical



Made in the U.S.A.

1610 West Jackson St. Unit 6
Painesville OH 44077
Office:440.350.9510 Fax:440.350.9520
service@faretec.com
www.FareTec.com

FareTec Inc.
Fracture and Rehabilitation Equipment