

RESPONDER GAUZE COMPRESSED COMBAT GAUZE



Responder Gauze takes care of moderate to severe wounds involving heavy blood loss. Made with high quality cotton for maximum absorption and stability. Responder Gauze is vacuum packed small for easy storage yet expands to cover large wound areas. It is the perfect companion to our compression bandages to stop severe blood loss and save lives.

This Gauze can be used to control hemorrhage in conjunction with a compression bandage, used as backing gauze for hemostatic dressings or for minor wound bandaging. Its low cube space requirement and rugged, durable vacuum sealed packaging make it easy to fit in both individual first aid kits and medic kits.

Special Features:

This 6-ply cotton gauze is conveniently z-folded which allows for a controlled application as a complete unit or in a continuous linear feed.

SPECIFICATIONS:

- MODEL # 2221038
- Package Dimensions: 6.3 x 5 x 1.27cm
- Package Weight: 1 oz.
- Case Contents: 260 units
- Shelf Life: 3 Years
- 5 mil poly bag, resists tears and holes
- Sterile, Z folded for easy application
- 10.4cm x 4.1 meters when opened

FareTec Inc Responder Products
1610 West Jackson St #6
Painesville OH 44077
USA
www.FareTec.com
www.ResponderProducts.com

