



RESPONDER REDI SPLINT FLEXIBLE ALUMINUM STAPLE

The Redi Splint bends into any simple curve increasing its strength for use on fractured limbs. Made of Volara foam and “dead soft” aluminum, this device is strong enough for pre- or post-cast care. Its light and compact enough to easily carry in a back pack or medic kit.



FEATURES:

- Cut to Size, customize the exact size if necessary using most scissors.
- Moldable aluminum, molds easily for various applications.
- Cleanable exterior, disinfect with standard cleaning supplier and/or bleach for reuse.
- Radiolucent material, x-rays and CT-scans can be taken without removing the splint.
- Weather resistant, functions across a wide range of temperatures.
- Lightweight, much lighter than plaster or fiberglass. Secure with a wrap or tape.
- Volara foam with closed cell, crosslinked technology

SPECIFICATIONS:

- CLOSED CELL FOAM THICKNESS: 8mm
- 4mm per side/8mm total.
- DEAD SOFT ALUMINUM THICKNESS: 1mm
- SIZES: 3 x 31”, 4” x 36”
- COLORS: BLACK, OLIVE DRAB
- WEIGHT: 3-4 ounces

