

RESPONDER PELVIC BINDER

EMERGENCY | RESCUE | EMS

High quality Pelvic Binder that is efficient and easy to use. Simplified strap system for fast application to apply necessary pressure in wrapping a patients fractured pelvis. Utilizes 1000D, flame retardant Cordura and MilSpec Velcro for grip strength and repeatability.

This Binder can be slid underneath the patient using various techniques. It can be easily cut to size or folded back on itself for the desired length and access gap. Intuitive velcro system is then placed across the patient and pulled tight with straps. Three easy steps.

Integrates with FareTec traction splints (as broken femurs are common with pelvic injuries).



FEATURES:

- Highly adjustable to fit various sized patients.
- Quick and easy intuitive strapping system.
- Vacuum packed small to fit in Kits.
- Integrates with the CT-6 & CT-7 Leg splints.
- Available in vacuum sealed package.
- Antimicrobial / Antibacterial
- X-Ray Translucent
- Made in USA
- Several sizes and color options
- FDA & CE/MDR CERTIFIED
- COLORS: BLACK, BLUE, CAMO, OTHERS

SPECIFICATIONS:

PART NUMBER: 8003000

SIZE: LENGTH - 58" 147cm
WIDTH - 8" 20cm

WEIGHT: 1 lb, .45 kg

MATERIAL: 1000D CORDURA
MilSpec Velcro, 2" NYLON

PACKED DIMS: 8 x 8 x "
[20 x 20 x 5 cm]

