



RESPONDER PELVIC BINDER

RESPONDER PELVIC BINDER

High quality Pelvic Binder that is efficient and easy to use. Uses simplified strap system for fast application to apply necessary pressure in wrapping a patients fractured pelvis. Moveable side pouch allows integration with FareTec leg traction splints. Utilizing 1000D, flame retardant Cordura and MilSpec Velcro for grip strength and repeatability.

This Binder can be slid underneath the patient using various techniques. It can be easily cut to size or folded back on itself for the desired length. Intuitive velcro system is then placed across the patient and pulled tight with straps. Three easy steps.

Femur fractures are also common with pelvic fracture injuries. This binder utilizes an adjustable pouch that integrates with FareTec leg traction splints.

Patents applied for*

SPECIFICATIONS:

PART NUMBER: 8003000

SIZE: LENGTH - 58" 147cm

WIDTH - 8" 20cm

WEIGHT: 1 lb, .45 kg

MATERIAL: 1000D CORDURA

MilSpec Velcro, 2" NYLON

PACKED DIMS: 8 x 8 x 3"

[20 x 20 x 7.5 cm]

PEDI SIZE AVAILABLE



FEATURES:

- Highly adjustable to fit various sized patients.
- Quick and easy intuitive strapping system.
- Vacuum packed small to fit in Kits.
- Integrates with the CT-6 & CT-7 Leg splints.
- Re-usable
- Antimicrobial / Antibacterial
- X-Ray Translucent
- Made in USA
- Several sizes and color options

