



FEATURED PRODUCTS

CT-7 LEG TRACTION SPLINT
ADULT / PEDIATRIC / BELOW KNEE APPLICATION

The CT-7 leg splint has been SIMPLIFIED; applies in 3 easy steps. It has a SHORTER footprint and no longer extends beyond the patient's foot. It has been made more rigid and indestructible with thicker carbon tubes. The traction mechanism is easy, fast and super precise by means of an acme threaded polycarbonate rod system. We are confident that this leg traction splint is the best choice for any medic.

NSN: 6515-01-706-4198. FDA, Made in USA.

*THE CT-6 IS ALSO STILL AVAILABLE



UPDATED CE & MDR



NEW RESPONDER PELVIC BINDER

Efficient and easy to use. Simplified strap system for fast application to apply necessary pressure in wrapping a patient's fractured pelvis. Allows integration with FareTec leg traction splints. Utilizing 1000D, flame retardant Cordura and MilSpec Velcro for grip strength and repeatability. Vacuum packaged.

- Highly adjustable to fit various sized patients.
- Antimicrobial / Antibacterial
- X-Ray Translucent
- Colors: BLACK, BLUE, CAMO
- FDA, CE & MDR, ISO
- Made in USA




RESPONDER COMPRESSED GAUZE

Upgraded gauze quality and packaging. Offered at a great low price with CE mark and necessary control documents. Red indicator tear notches now included. Includes instruction packet with 26 languages.

RESPONDER COMPRESSION BANDAGES

Upgrade packaging on all bandages. Extra durable and double layered to protect the product and keep it vacuum sealed in the harshest environments. Features unique Manufacturer, expiration and LOT#s.



4" & 6" Sizes +ABD and FLAT fold



FARETEC INC

RESPONDER PRODUCTS

PORTABLE SURGERY TABLE (PST) PATIENT TRANSPORT AND OPERATION PLATFORM

Field tested and outfitted with the most advanced materials available, this unit can serve a multitude of purposes. The PST quickly assembles into a stable patient aid platform for medics and doctors. The device is versatile and easy to transport for fast changing scenarios.


Standard Accessories:

- IV Pole [28-52in./70-132cm] Height Adjustable.
- Mayo Tray [13x9in./33x23cm] Height and Rotation Adjustable.
- Armboard, Padded [3/8" thick black phenolic, 4.5 lbs] Radiolucent,
- LED Exam Light, Battery operated. Fully Adjustable. 24 hr life.
- Storage Catcher, Convenient, under patient nylon fabric



(IFAK) INDIVIDUAL FIRST AID KIT & BACKPACKS

Top quality bags combined with the best life saving products available, the Responder IFAK has many advantages. Each unit is packaged small enough to be carried anywhere, rugged enough for even the most extreme environments, sterile and ready for immediate use with a minimum amount of training. We are an authorized distributor of the following items:
CATGEN7 TOURNIQUETS, CELOX HEMOSTAT. FOX & HALO CHEST SEALS, TYTEK 14GA NEEDLES



*The contents of each IFAK can be altered to fit customer needs; we build to your specifications.

LIFEBOX50 TRANSPORT CONTAINERS CLIMATE CONTROLLED

LifeBox50 is intended for the protection, transport and temporary isolation of temperature-sensitive content such as vaccines, medications, cold therapy items, and blood. Rugged Pelican Hardware combined with R50 insulation rating.

TWO NEW SIZES AVAILABLE:

- LB50 BRIEFCASE / Part # 6001023 - INNER DIMENSIONS 8.25 x 3.25 x 13.5in
- SMALL LB50 NSN#: 8145 01 596 4629
- LB50 XSMALL / Part # 60001022 - INNER DIMENSIONS 6.75 x 3.5 x 10.5in



Each unit is mounted with a digital thermometer and includes reusable cold gel packs.