



THE C-A-TOURNIQUET®

**1ST CHOICE
OF THE U.S. MILITARY**

The **COMBAT APPLICATION TOURNIQUET®** is a small one-handed tourniquet that completely occludes arterial blood flow in an extremity.

The **C-A-T®** uses a Self-Adhering Band and a Friction Adaptor Buckle to fit a wide range of extremities combined with a one-handed Windlass System™.

The Windlass™ uses a free moving internal band to provide true circumferential pressure to an extremity.

The Windlass™ is then locked in place; this requires only one hand with the Windlass Clip™.

**#1 LIFE
SAVING
TOURNIQUET**



**COMPOSITE
RESOURCES**



C-A-TOURNIQUET®

The **C-A-TOURNIQUET®** is proven to be 100% effective in occluding blood flow in both upper & lower extremities by the U.S. Army's Institute of Surgical Research.

The **C-A-T®** also has a Hook-and-Loop Windlass Strap™ for further securing of the Windlass™ during patient transport. The latest generation of the **C-A-TOURNIQUET®** employs unique features:

- True One-Handed Application
- Reinforced Windlass Clip™
- Quickly controls bleeding
- Highly Visibly Time Strap
- Enhanced Windlass Rod™
- Strong & Lightweight