

## QD-4 HARE ADULT & PEDIATRIC TRACTION SPLINT

Used by emergency services across the country, the QD line of splint kits has several unique and proprietary features. The “Hex Lock” system allows the splint to be quickly adjusted to any length. Once it has been sized appropriately, the intuitive locking mechanism is engaged by turning the Hex Lock to either side. It’s fast, easy and reliable. This system helps eliminate problems with dirt or corrosion which cause other splints to malfunction. It has become an industry standard for patients requiring Hare style traction splinting.

### FEATURES:

- Compact, Strong, Light.
- Aluminum or Molded Ratchet.
- Integral Strap System.
- Field tested.
- QD-4 fits Pedi to Adult.
- Intuitive Adjusting System.
- Minimal patient movement.
- QD-3 Pediatric traction splint for states that require separate unit.



### SPECIFICATIONS:

- QD-4 Length Extended 54" (137cm)  
Length Closed 32" (81cm)
- QD-3 Length Extended 45" (114cm)  
Length Closed 29" (73cm)
- Weight: 10.4 lbs (4.7 kilos)



The QD-4 splint will adjust to accommodate pediatrics and adults. All units come as a complete kit with straps, ankle hitch and carrying case included. In addition, all QD's are available with either our anodized aluminum ratchet or our injection molded ratchet.



| Part Number | Description           | Ratchet  |
|-------------|-----------------------|----------|
| 1126524     | QD-4 Adult Splint     | Aluminum |
| 1126522     | QD-3 Pediatric Splint | Aluminum |
| 1126514     | QD-4 Adult Splint     | Molded   |
| 1126512     | QD-3 Pediatric Splint | Molded   |
| 1126726     | QD 3&4 Combination    | Aluminum |
| 1126716     | QD 3&4 Combination    | Molded   |