





## RESPONDER IFAK





ITEM NAME	DESCRIPTION
<p>IFAK BAG</p> 	<p>Durable design, construction and stitching. Made of rugged 1000 Denier. Double Zipper top to side opening. Internal elastic attachment points and Elastic straps to prevent over-opening. MOLLE, ALICE and belt combatability. Color: KSA military camo.</p> <p>Dimensions: 21 x 14 x 6cm Factory of gabardine high quality prevents the entry of air and is anti-flammable. Conforms to international standards. Internal pockets to hold multiple items with elastic; even when upside down. Outside pocket for quick access to essential items.</p>
<p>CAT TOURNIQUET</p> 	<p>The patented C-A-T® is a true one-handed tourniquet proven to be 100% effective by the U.S. Army's Institute of Surgical Research. Tests proved that the C-A-T® completely occluded blood flow of an extremity in the event of a traumatic wound with significant hemorrhage. The C-A-T® utilizes a durable windlass system with a patented free-moving internal band providing true circumferential pressure to the extremity. Once adequately tightened, bleeding will cease and the windlass is locked into place. A hook and loop windlass retention strap is then applied, securing the windlass to maintain pressure during casualty evacuation. The C-A-T®'s unique dual securing system avoids the use of screws and clips which can become difficult to operate under survival stress or where fine motor skills are compromised.</p>





## RESPONDER IFAK



ITEM NAME	DESCRIPTION
<p>RESPONDER COMPRESSED GAUZE</p> 	<p>Responder Compressed Gauze is an ultra-compact roll of sterile, high quality 100% cotton gauze that utilizes a unique crinkle weave that has excellent fluid absorption. It is a 6-ply, 4.5" x 4.1 yd. gauze with finished edges to reduce loose ends or lint. The gauze is indicated for the treatment of bleeding and the dressing and bandaging of wounds. Responder Gauze features a proprietary packaging process that ensures the most consistent, low-cube package on the market. The package is a rugged and durable vacuum sealed container with easy-to-open tear notches. The ultra-compact storage size expands once opened to cover larger wound areas.</p>
<p>RESPONDER COMPRESSION BANDAGE</p> 	<p>The Responder Bandage is designed for quick and easy application by caregivers to provide effective, multi-functional treatment when time is of the importance.</p> <p>The Responder Bandage's sterile, non-adherent pad applies pressure to any site, is easily wrapped and secured. The Responder Bandage will restrict blood flow to help in the treatment of Traumatic Bleeds.</p> <p>The Responder Bandage is available in 2 sizes, 4" and 6", specifications are:  4" bandage  Pad wt: three layer, 280g+280g+110g/m2  Pad size: 10x18cm  Bandage wt: 80-85g/m2  Width of bandage: 10cm  Natural length: 1.5m  Elas_c length: 4.5m  Double bag rolling packing, inner bag vacuum packing.</p>





## RESPONDER IFAK



ITEM NAME	DESCRIPTION
<p>ADHESIVE TAPE SURGICAL ROLL 2 INCH</p>	<p>A great item that can be used for multiple tasks in order to secure and combine items to the patient and each other. WHITE</p>
<p>NITRILE GLOVES, LG</p>	<p>Non-latex, powder-free nitrile barrier gloves are an excellent option for the many people for whom latex causes severe allergic reactions and also provides increased protection against viruses. PAIR of 2, STERILE</p>
<p>TACTICAL COMBAT CARD ARABIC</p>	<p>Quickly references injuries a patient has previously received and Identifies the patient. This card allows next responder to know what medications have been given along with any other pertinent information.</p>
<p>F28 NASAL PHARYNGEAL &amp; LUBE RESPONDER NP</p>	<p>It is used to maintain an open airway, by inserting the tube into the nasal passageway. When a patient becomes unconscious, the muscles in the jaw commonly relax and can allow the tongue to slide back and obstruct the airway. The purpose of the flared end is to prevent the device from becoming lost inside the patient's nose.</p> <ul style="list-style-type: none"> <li>- Maximum kink-resistance helps to ensure a patent airway</li> <li>- Soft PVC, nontoxic, flexible design</li> <li>- Smoothly finished and rounded edges, less nasal trauma, maximize patient comfort</li> <li>- Smooth surface facilitates insertion</li> <li>- Size identified on the tube</li> <li>- Tube with flared end preventing the device from becoming lost inside the patient's nose</li> <li>- Latex Free</li> </ul>









## RESPONDER IFAK



ITEM NAME	DESCRIPTION
<p>TRAUMA SHEERS</p> 	<p>These paramedic shears have a serrated lower edge which can cut through the toughest materials. Their soft-grip, non-slip contoured handles provide maximum control, even when wet. Their clothing lift lip and rounded tip provide the ultimate in patient safety. These black shears are ideal for use in tactical settings and are made to meet the challenges of casualty care with precision and dependability.</p>
<p>QUIKLOT HEMOSTATIC GAUZE</p> 	<p>QuikClot Combat Gauze Z-Fold is a soft, white, sterile, nonwoven gauze impregnated with kaolin, an inert mineral that does not contain animal or human proteins or botanicals. QuikClot Combat Gauze Z-Fold is an improvement on the original QuikClot Combat Gauze. As a result feedback from the US Military, Z-Medica improved the QuikClot Combat Gauze product by making it z-folded, vacuum packed, and with an x-ray detectable strip. The QuikClot Combat Gauze Z-Fold is even easier to apply since the gauze is in a z-fold format, which makes packing wounds easier and eliminates the risk of the roll falling away and getting contaminated. The vacuum packed feature of the QuikClot Combat Gauze Z-Fold makes it much more compact, saving room on your body. Finally, the x-ray detectable strip that is incorporated in the QuikClot Combat Gauze Z-Fold ensures that the gauze will be found if the patient is taken to the operating room.</p>